

How will we determine if Calcineurin binds a specific NFAT fragment?

You have:

- Cell lysate from bacteria expressing GST or GST-NFAT fragment
- Concentrated glutathione
- Glutathione-agarose columns
- Purified Calcineurin protein
- SDS buffer
- A buffer containing calcium

Design an experiment to test if Calcineurin binds your NFAT fragment

How will we determine if Calcineurin binds a specific NFAT fragment?

You have:

- 1 - Cell lysate from bacteria expressing GST or GST-NFAT fragment
- 4 - Concentrated glutathione
- 2 - Glutathione-agarose columns ←
- 3 - Purified Calcineurin protein
- 5 - SDS buffer → SDS-PAGE → western blot
- 3 - A buffer containing calcium

Design an experiment to test if Calcineurin binds your NFAT fragment

Design an experiment to test if Calcineurin binds your NFAT fragment

